



Date d'édition : 18.05.2024

Ref : 8-1004196-100-10-0

P 4.196 Lampe à vapeur de sodium à haute pression

P 4.196 Experimental Panel: High Pressure Sodium Vapour Lamp 100 W Ballast device and electronic superimposed pulse igniter for a safe lamp start. An automatic cut-out safely disconnects the lamp after 82 seconds at end of its service life or in the event of a defect. - Igniter voltage: 4.0 ... 5.0 kV - Compensation condenser: 12 μ F / 400 V - Lamp socket with ballast and starter: E 40, 230 V AC, 50 ... 60 Hz, 100 W. 1.2 A Dimensions in mm: 320 x 297 x 255 (W x H x D) Mass: 3.36 kg