



Date d'édition : 19.05.2026

**Ref : BEDK-004**

**Bacterial Plasmid DNA**

Analysis of bacterial plasmid DNA in gel electrophoresis. Here, for example, uncut (circular) plasmid DNA can be compared to plasmid DNA that has already been enzymatically cut (linearised). Sufficient for 10 individual tests.

Matériel livré :

DNA uncut  
Cut DNA with EcoRI, HinfI and HaeIII  
TAE electrophoresis buffer  
Agarose  
DNA staining concentrate  
Test description in German

En option:

Required accessories for pouring the DNA gel: Erlenmeyer flask, microwave, lab scale and dist. water. Storage currently at 4°C (-18°C for long-term storage).

Catégories / Arborescence

Sciences > Biologie > Produits > Genetics and evolution